

Medical Antipyretic Patch


Description:

An antipyretic patch, also known as a fever patch, fever cooling gel patch or cooling patch is a medical product designed to provide relief from fever and associated discomfort. Medical antipyretic patches are topical applications designed to provide localized cooling to the skin. They are commonly used for pain relief, reducing inflammation, and managing fever.



Features:

Convenient, Portable, Disposable
Place Firmly, Not Easy to Fall



Rapid Onset:

Cooling patches provide an almost immediate cooling sensation upon application, which can help alleviate discomfort quickly.

High Water Content:

The hydrogel layer of medical antipyretic patches contains a high percentage of water, which evaporates to produce a cooling effect.

Flexible and Soft:

Hydrogels are soft and flexible, conforming to the contours of the body for comfortable wear.

Skin-Friendly:

Medical-grade adhesives are used to ensure the patch stays in place without causing irritation or discomfort.

Hypoallergenic:

Cooling patches are designed to minimize the risk of allergic reactions, making them suitable for sensitive skin.

Ease of Use:

Easy to apply and remove antipyretic patches, making them convenient for use at home, work, or on-the-go.

Long-Lasting:

Our cooling patches are designed to provide relief for 8 hours, reducing the need for frequent reapplication.


Portable:

Medical antipyretic patches are often individually packaged, making them easy to carry in a bag or pocket for use whenever needed.

Specification

Product Name:	Medical Antipyretic Patch/ Cooling Gel Patch/ Fever Patch / Cooling Pad
Patch size:	5*12cm, 5*15cm, 4*8cm, 4*10 cm or customized
Gel color:	Blue or customized
Feature	Physical cooling
Material:	Non-woven fabric back lining, gel layer, PE protective film
Application	Apply to adult and children's body parts to cool down

Material



The medical antipyretic patch is made of non-woven fabric back lining, gel layer and PE protective film. The Main Ingredients include hydrogel, purified water, peppermint, pigment. The product does not include medical ingredients related with pharmacology, immunology or metabolism.

Application:



Fever



Feel Sleepy




Hot Day

Instruction:



Step 1:

Tear off the transparent film of the medical antipyretic patch.



Step 2:

Apply the medical antipyretic patch to the area that needs cooling.

